

## Product datasheet for **TL312084**

### CD11c (ITGAX) Human shRNA Plasmid Kit (Locus ID 3687)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	CD11c (ITGAX) Human shRNA Plasmid Kit (Locus ID 3687)
Locus ID:	3687
Synonyms:	CD11C; SLEB6
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	ITGAX - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 3687). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_000887</a> , <a href="#">NM_001286375</a> , <a href="#">NM_000887.1</a> , <a href="#">NM_000887.2</a> , <a href="#">NM_000887.3</a> , <a href="#">NM_000887.4</a> , <a href="#">NM_001286375.1</a> , <a href="#">BC037577</a> , <a href="#">BC038237</a> , <a href="#">BM149283</a> , <a href="#">NM_000887.5</a> , <a href="#">NM_001286375.2</a>
UniProt ID:	<a href="#">P20702</a>
Summary:	This gene encodes the integrin alpha X chain protein. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. This protein combines with the beta 2 chain (ITGB2) to form a leukocyte-specific integrin referred to as inactivated-C3b (iC3b) receptor 4 (CR4). The alpha X beta 2 complex seems to overlap the properties of the alpha M beta 2 integrin in the adherence of neutrophils and monocytes to stimulated endothelium cells, and in the phagocytosis of complement coated particles. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).