

## Product datasheet for **TL312114**

### IREB2 Human shRNA Plasmid Kit (Locus ID 3658)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	IREB2 Human shRNA Plasmid Kit (Locus ID 3658)
Locus ID:	3658
Synonyms:	ACO3; IRE-BP 2; IRE-BP2; IRP2; IRP2AD; NDCAMA
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	IREB2 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 3658). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001320941</a> , <a href="#">NM_001320942</a> , <a href="#">NM_001320943</a> , <a href="#">NM_004136</a> , <a href="#">NM_001354994</a> , <a href="#">NM_004136.1</a> , <a href="#">NM_004136.2</a> , <a href="#">NM_004136.3</a> , <a href="#">BC017880</a> , <a href="#">BC117481</a> , <a href="#">BC117483</a> , <a href="#">BC143948</a> , <a href="#">NM_004136.4</a>
UniProt ID:	<a href="#">P48200</a>
Summary:	The protein encoded by this gene is an RNA-binding protein that acts to regulate iron levels in the cells by regulating the translation and stability of mRNAs that affect iron homeostasis under conditions when iron is depleted. When iron levels are low, this protein binds to iron-responsive elements (IRES), stem-loop structures located either in the 5' or 3' UTRs. Binding to the 5' UTR represses translation, while binding to the 3' UTR inhibits mRNA degradation. When iron is found in the cell, this protein is degraded in a F-box and leucine rich repeat protein 5-dependent manner. Variants in this gene have been associated with lung cancer and chronic obstructive pulmonary disease (COPD). Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Aug 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).