

Product datasheet for **TL312123**

IQ motif containing B1 (IQCB1) Human shRNA Plasmid Kit (Locus ID 9657)

Product data:

Product Type:	shRNA Plasmids
Product Name:	IQ motif containing B1 (IQCB1) Human shRNA Plasmid Kit (Locus ID 9657)
Locus ID:	9657
Synonyms:	NPHP5; PIQ; SLSN5
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	<p>IQCB1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 9657). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.</p>
RefSeq:	NM_001023570 , NM_001023571 , NM_001319107 , NM_014642 , NR_134968 , NM_001023571.1 , NM_001023571.2 , NM_001023571.3 , NM_001023570.1 , NM_001023570.2 , NM_001023570.3 , BC005806 , BC005806.2 , BM682018 , NM_001023570.4
UniProt ID:	Q15051
Summary:	<p>This gene encodes a nephrocystin protein that interacts with calmodulin and the retinitis pigmentosa GTPase regulator protein. The encoded protein has a central coiled-coil region and two calmodulin-binding IQ domains. It is localized to the primary cilia of renal epithelial cells and connecting cilia of photoreceptor cells. The protein is thought to play a role in ciliary function. Defects in this gene result in Senior-Loken syndrome type 5. Alternative splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 6. [provided by RefSeq, Jan 2016]</p>
shRNA Design:	<p>These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our custom shRNA service.</p>


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).