

## Product datasheet for **TL312185**

### IL18R Beta (IL18RAP) Human shRNA Plasmid Kit (Locus ID 8807)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	IL18R Beta (IL18RAP) Human shRNA Plasmid Kit (Locus ID 8807)
Locus ID:	8807
Synonyms:	ACPL; CD218b; CDw218b; IL-1R-7; IL-1R7; IL-1RAcPL; IL-18R-beta; IL-18RAcP; IL-18Rbeta; IL18RB
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	IL18RAP - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 8807). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_003853</a> , <a href="#">NM_003853.1</a> , <a href="#">NM_003853.2</a> , <a href="#">NM_003853.3</a> , <a href="#">BC106765</a> , <a href="#">BC069630</a> , <a href="#">BC106764</a> , <a href="#">BC137474</a> , <a href="#">BC137475</a>
UniProt ID:	<a href="#">O95256</a>
Summary:	The protein encoded by this gene is an accessory subunit of the heterodimeric receptor for interleukin 18 (IL18), a proinflammatory cytokine involved in inducing cell-mediated immunity. This protein enhances the IL18-binding activity of the IL18 receptor and plays a role in signaling by IL18. Mutations in this gene are associated with Crohn's disease and inflammatory bowel disease, and susceptibility to celiac disease and leprosy. Alternatively spliced transcript variants of this gene have been described, but their full-length nature is not known. [provided by RefSeq, Feb 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).