

## Product datasheet for **TL312186**

### IL18R1 Human shRNA Plasmid Kit (Locus ID 8809)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	IL18R1 Human shRNA Plasmid Kit (Locus ID 8809)
Locus ID:	8809
Synonyms:	CD218a; CDw218a; IL-1Rrp; IL-18R-alpha; IL-18Ralpha; IL1RRP; IL18RA; IL18Ralpha2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	IL18R1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 8809). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001282399</a> , <a href="#">NM_003855</a> , <a href="#">NM_003855.1</a> , <a href="#">NM_003855.2</a> , <a href="#">NM_003855.3</a> , <a href="#">NM_001282399.1</a> , <a href="#">BC069575</a> , <a href="#">BC069575.1</a> , <a href="#">BC093975</a> , <a href="#">BC093977</a> , <a href="#">BC143403</a> , <a href="#">BC143405</a> , <a href="#">NM_003855.5</a> , <a href="#">NM_001282399.2</a>
UniProt ID:	<a href="#">Q13478</a>
Summary:	The protein encoded by this gene is a cytokine receptor that belongs to the interleukin 1 receptor family. This receptor specifically binds interleukin 18 (IL18), and is essential for IL18 mediated signal transduction. IFN-alpha and IL12 are reported to induce the expression of this receptor in NK and T cells. This gene along with four other members of the interleukin 1 receptor family, including IL1R2, IL1R1, ILRL2 (IL-1Rrp2), and IL1RL1 (T1/ST2), form a gene cluster on chromosome 2q. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).