

## Product datasheet for **TL312212**

### CD101 Human shRNA Plasmid Kit (Locus ID 9398)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	CD101 Human shRNA Plasmid Kit (Locus ID 9398)
Locus ID:	9398
Synonyms:	EWI-101; IGSF2; V7
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	CD101 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 9398). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001256106</a> , <a href="#">NM_001256109</a> , <a href="#">NM_001256111</a> , <a href="#">NM_004258</a> , <a href="#">NM_004258.1</a> , <a href="#">NM_004258.2</a> , <a href="#">NM_004258.3</a> , <a href="#">NM_004258.4</a> , <a href="#">NM_004258.5</a> , <a href="#">NM_001256106.1</a> , <a href="#">NM_001256106.2</a> , <a href="#">NM_001256109.1</a> , <a href="#">NM_001256109.2</a> , <a href="#">NM_001256111.1</a> , <a href="#">NM_001256111.2</a> , <a href="#">BC130327</a> , <a href="#">BC136484</a> , <a href="#">BC171768</a>
UniProt ID:	<a href="#">Q93033</a>
Summary:	Plays a role as inhibitor of T-cells proliferation induced by CD3. Inhibits expression of IL2RA on activated T-cells and secretion of IL2. Inhibits tyrosine kinases that are required for IL2 production and cellular proliferation. Inhibits phospholipase C-gamma-1/PLCG1 phosphorylation and subsequent CD3-induced changes in intracellular free calcium. Prevents nuclear translocation of nuclear factor of activated T-cell to the nucleus. Plays a role in the inhibition of T-cell proliferation via IL10 secretion by cutaneous dendritic cells. May be a marker of CD4(+) CD56(+) leukemic tumor cells.[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).