

Product datasheet for **TL312222V**

IGF2BP2 Human shRNA Lentiviral Particle (Locus ID 10644)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	IGF2BP2 Human shRNA Lentiviral Particle (Locus ID 10644)
Locus ID:	10644
Synonyms:	IMP-2; IMP2; VICKZ2
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	IGF2BP2 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	<u>NM_001007225</u> , <u>NM_001291869</u> , <u>NM_001291872</u> , <u>NM_001291873</u> , <u>NM_001291874</u> , <u>NM_001291875</u> , <u>NM_006548</u> , <u>NR_138486</u> , <u>NM_006548.2</u> , <u>NM_006548.3</u> , <u>NM_006548.4</u> , <u>NM_006548.5</u> , <u>NM_001007225.1</u> , <u>NM_001007225.2</u> , <u>NM_001291875.1</u> , <u>NM_001291875.2</u> , <u>NM_001291874.1</u> , <u>NM_001291874.2</u> , <u>NM_001291873.1</u> , <u>NM_001291873.2</u> , <u>NM_001291872.1</u> , <u>NM_001291872.2</u> , <u>NM_001291869.1</u> , <u>NM_001291869.2</u> , <u>BC021290</u> , <u>BC021290.2</u> , <u>NM_001291875.3</u> , <u>NM_001007225.3</u> , <u>NM_001291874.3</u> , <u>NM_001291869.3</u> , <u>NM_001291872.3</u>
UniProt ID:	<u>Q9Y6M1</u>
Summary:	This gene encodes a protein that binds the 5' UTR of insulin-like growth factor 2 (IGF2) mRNA and regulates its translation. It plays an important role in metabolism and variation in this gene is associated with susceptibility to diabetes. Alternative splicing and promoter usage results in multiple transcript variants. Related pseudogenes are found on several chromosomes. [provided by RefSeq, Sep 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).