

Product datasheet for **TL312338**

HRH1 Human shRNA Plasmid Kit (Locus ID 3269)

Product data:

Product Type:	shRNA Plasmids
Product Name:	HRH1 Human shRNA Plasmid Kit (Locus ID 3269)
Locus ID:	3269
Synonyms:	H1-R; H1R; HH1R; hisH1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	HRH1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 3269). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_000861 , NM_001098211 , NM_001098212 , NM_001098213 , NM_000861.1 , NM_000861.2 , NM_000861.3 , NM_001098213.1 , NM_001098212.1 , NM_001098211.1 , BC060802 , BC060802.1 , BM905371 , NM_001098211.2 , NM_001098213.2 , NM_001098212.2
UniProt ID:	P35367
Summary:	Histamine is a ubiquitous messenger molecule released from mast cells, enterochromaffin-like cells, and neurons. Its various actions are mediated by histamine receptors H1, H2, H3 and H4. The protein encoded by this gene is an integral membrane protein and belongs to the G protein-coupled receptor superfamily. It mediates the contraction of smooth muscles, the increase in capillary permeability due to contraction of terminal venules, the release of catecholamine from adrenal medulla, and neurotransmission in the central nervous system. It has been associated with multiple processes, including memory and learning, circadian rhythm, and thermoregulation. It is also known to contribute to the pathophysiology of allergic diseases such as atopic dermatitis, asthma, anaphylaxis and allergic rhinitis. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Jan 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL312338A virus into HEK293 cells. TL312338A virus was prepared using lenti-shRNA TL312338A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL312338B virus into HEK293 cells. TL312338B virus was prepared using lenti-shRNA TL312338B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL312338C] virus into HEK293 cells. [TL312338C] virus was prepared using lenti-shRNA [TL312338C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL312338D] virus into HEK293 cells. [TL312338D] virus was prepared using lenti-shRNA [TL312338D] and [TR30037] packaging kit.