

Product datasheet for **TL312404V**

HLA-DRB1 Human shRNA Lentiviral Particle (Locus ID 3123)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	HLA-DRB1 Human shRNA Lentiviral Particle (Locus ID 3123)
Locus ID:	3123
Synonyms:	DRB1; HLA-DR1B; HLA-DRB; SS1
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	HLA - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	BC007920 , BC008403 , BC024269 , BC033827 , NM_001243965 , NM_002124 , NM_001359193 , NM_001359194 , NM_002124.1 , NM_002124.2 , NM_002124.3 , NM_001243965.1 , BC007920.2 , BC024269.1 , BC033827.1 , BC018832 , BC018834 , BC018835 , BC031023 , BC108922 , BM671866
UniProt ID:	P04229
Summary:	HLA-DRB1 belongs to the HLA class II beta chain paralogs. The class II molecule is a heterodimer consisting of an alpha (DRA) and a beta chain (DRB), both anchored in the membrane. It plays a central role in the immune system by presenting peptides derived from extracellular proteins. Class II molecules are expressed in antigen presenting cells. The beta chain is approximately 26-28 kDa. It is encoded by 6 exons. Exon one encodes the leader peptide; exons 2 and 3 encode the two extracellular domains; exon 4 encodes the transmembrane domain; and exon 5 encodes the cytoplasmic tail. Within the DR molecule the beta chain contains all the polymorphisms specifying the peptide binding specificities. Hundreds of DRB1 alleles have been described and some alleles have increased frequencies associated with certain diseases or conditions. For example, DRB1*1302 has been related to acute and chronic hepatitis B virus persistence. There are multiple pseudogenes of this gene. [provided by RefSeq, Jul 2020]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL312404A virus into HEK293 cells. TL312404A virus was prepared using lenti-shRNA TL312404A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL312404B virus into HEK293 cells. TL312404B virus was prepared using lenti-shRNA TL312404B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL312404C] virus into HEK293 cells. [TL312404C] virus was prepared using lenti-shRNA [TL312404C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL312404D] virus into HEK293 cells. [TL312404D] virus was prepared using lenti-shRNA [TL312404D] and [TR30037] packaging kit.