

## Product datasheet for **TL312414**

### ZBTB48 Human shRNA Plasmid Kit (Locus ID 3104)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	ZBTB48 Human shRNA Plasmid Kit (Locus ID 3104)
Locus ID:	3104
Synonyms:	HKR3; pp9964; TZAP; ZNF855
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	ZBTB48 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 3104). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001278647</a> , <a href="#">NM_001278648</a> , <a href="#">NM_005341</a> , <a href="#">NM_005341.1</a> , <a href="#">NM_005341.2</a> , <a href="#">NM_005341.3</a> , <a href="#">NM_001278647.1</a> , <a href="#">NM_001278648.1</a> , <a href="#">BC013573</a> , <a href="#">NM_001278648.2</a> , <a href="#">NM_001278647.2</a> , <a href="#">NM_005341.4</a>
UniProt ID:	<a href="#">P10074</a>
Summary:	Telomere-binding protein that acts as a regulator of telomere length (PubMed:28500257, PubMed:28082411). Directly binds the telomeric double-stranded 5'-TTAGGG-3' repeat (PubMed:28500257, PubMed:28082411). Preferentially binds to telomeres that have a low concentration of shelterin complex and acts as a regulator of telomere length by initiating telomere trimming, a process that prevents the accumulation of aberrantly long telomeres (PubMed:28082411). Also acts as a transcription regulator that binds to promoter regions (PubMed:7969177, PubMed:24382891, PubMed:28500257). Regulates expression of a small subset of genes, including MTFP1 (PubMed:28500257). Regulates expression the J and/or S elements in MHC II promoter (PubMed:7969177). Acts as a negative regulator of cell proliferation by specifically activating expression of ARF, a tumor suppressor isoform of CDKN2A (PubMed:24382891).[UniProtKB/Swiss-Prot Function]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).