

## Product datasheet for **TL312704**

### Golgin subfamily A member 4 (GOLGA4) Human shRNA Plasmid Kit (Locus ID 2803)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Golgin subfamily A member 4 (GOLGA4) Human shRNA Plasmid Kit (Locus ID 2803)
Locus ID:	2803
Synonyms:	CRPF46; GCP2; GOLG; MU-RMS-40.18; p230
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	GOLGA4 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 2803). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001172713</a> , <a href="#">NM_002078</a> , <a href="#">NM_002078.1</a> , <a href="#">NM_002078.2</a> , <a href="#">NM_002078.3</a> , <a href="#">NM_002078.4</a> , <a href="#">NM_001172713.1</a> , <a href="#">BC023648</a> , <a href="#">BC045543</a> , <a href="#">BC050444</a> , <a href="#">BC172405</a> , <a href="#">NM_002078.5</a> , <a href="#">NM_001172713.2</a>
UniProt ID:	<a href="#">Q13439</a>
Summary:	The Golgi apparatus, which participates in glycosylation and transport of proteins and lipids in the secretory pathway, consists of a series of stacked cisternae (flattened membrane sacs). Interactions between the Golgi and microtubules are thought to be important for the reorganization of the Golgi after it fragments during mitosis. This gene encodes one of the golgins, a family of proteins localized to the Golgi. This protein has been postulated to play a role in Rab6-regulated membrane-tethering events in the Golgi apparatus. Alternatively spliced transcript variants encoding different isoforms have been identified in this gene. [provided by RefSeq, Feb 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).