

Product datasheet for **TL312779**

Growth hormone receptor (GHR) Human shRNA Plasmid Kit (Locus ID 2690)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Growth hormone receptor (GHR) Human shRNA Plasmid Kit (Locus ID 2690)
Locus ID:	2690
Synonyms:	GHBP; GHIP
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	GHR - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 2690). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_000163 , NM_001242399 , NM_001242400 , NM_001242401 , NM_001242402 , NM_001242403 , NM_001242404 , NM_001242405 , NM_001242406 , NM_001242460 , NM_001242461 , NM_001242462 , NM_000163.1 , NM_000163.2 , NM_000163.3 , NM_000163.4 , NM_001242462.1 , NM_001242461.1 , NM_001242460.1 , NM_001242400.1 , NM_001242400.2 , NM_001242401.1 , NM_001242401.2 , NM_001242401.3 , NM_001242402.1 , NM_001242402.2 , NM_001242403.1 , NM_001242403.2 , NM_001242404.1 , NM_001242404.2 , NM_001242405.1 , NM_001242405.2 , NM_001242406.1 , NM_001242406.2 , NM_001242399.1 , NM_001242399.2 , BC136496 , BM719293 , NM_001242401.4 , NM_000163.5
UniProt ID:	P10912



[View online »](#)

- Summary:** This gene encodes a member of the type I cytokine receptor family, which is a transmembrane receptor for growth hormone. Binding of growth hormone to the receptor leads to receptor dimerization and the activation of an intra- and intercellular signal transduction pathway leading to growth. Mutations in this gene have been associated with Laron syndrome, also known as the growth hormone insensitivity syndrome (GHIS), a disorder characterized by short stature. In humans and rabbits, but not rodents, growth hormone binding protein (GHBP) is generated by proteolytic cleavage of the extracellular ligand-binding domain from the mature growth hormone receptor protein. Multiple alternatively spliced transcript variants have been found for this gene.[provided by RefSeq, Jun 2011]
- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).