

Product datasheet for **TL312854**

GAL3ST1 Human shRNA Plasmid Kit (Locus ID 9514)

Product data:

Product Type:	shRNA Plasmids
Product Name:	GAL3ST1 Human shRNA Plasmid Kit (Locus ID 9514)
Locus ID:	9514
Synonyms:	CST
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	GAL3ST1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 9514). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001318103 , NM_001318104 , NM_001318105 , NM_001318106 , NM_001318107 , NM_001318108 , NM_001318109 , NM_001318110 , NM_001318111 , NM_001318112 , NM_001318113 , NM_001318114 , NM_001318115 , NM_001318116 , NM_004861 , NM_004861.2 , BC019077 , BC019077.2 , BC014649 , BM923857 , NM_004861.3
UniProt ID:	Q99999
Summary:	Sulfonation, an important step in the metabolism of many drugs, xenobiotics, hormones, and neurotransmitters, is catalyzed by sulfotransferases. This gene encodes galactosylceramide sulfotransferase, which catalyzes the sulfation of membrane glycolipids including the final step in the synthesis of sulfatide, a major lipid component of the myelin sheath. This gene exhibits elevated expression in ovarian epithelial carcinoma and the encoded enzyme exhibits elevated activity in renal cell carcinoma. Mutations in this gene may be associated with reduced insulin resistance. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2015]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).