

Product datasheet for **TL312893**

Fuc-TIX (FUT9) Human shRNA Plasmid Kit (Locus ID 10690)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Fuc-TIX (FUT9) Human shRNA Plasmid Kit (Locus ID 10690)
Locus ID:	10690
Synonyms:	Fuc-TIX
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	FUT9 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 10690). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_006581 , NM_006581.1 , NM_006581.2 , NM_006581.3 , BC036101 , BC001879 , NM_006581.4
UniProt ID:	Q9Y231
Summary:	The protein encoded by this gene belongs to the glycosyltransferase family. It is localized to the golgi, and catalyzes the last step in the biosynthesis of Lewis X (LeX) antigen, the addition of a fucose to precursor polysaccharides. This protein is one of the few fucosyltransferases that synthesizes the LeX oligosaccharide (CD15) expressed in the organ buds progressing in mesenchyma during embryogenesis. It is also responsible for the expression of CD15 in mature granulocytes. A common haplotype of this gene has also been associated with susceptibility to placental malaria infection. [provided by RefSeq, Nov 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).