

## Product datasheet for **TL312965**

### FLRT2 Human shRNA Plasmid Kit (Locus ID 23768)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	FLRT2 Human shRNA Plasmid Kit (Locus ID 23768)
Locus ID:	23768
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	FLRT2 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 23768). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001346143</a> , <a href="#">NM_001346144</a> , <a href="#">NM_001346145</a> , <a href="#">NM_001346146</a> , <a href="#">NM_013231</a> , <a href="#">NR_144385</a> , <a href="#">NR_144386</a> , <a href="#">NR_144387</a> , <a href="#">NM_013231.1</a> , <a href="#">NM_013231.2</a> , <a href="#">NM_013231.3</a> , <a href="#">NM_013231.4</a> , <a href="#">NM_013231.5</a> , <a href="#">BC126249</a> , <a href="#">BC130290</a> , <a href="#">BC143936</a>
UniProt ID:	<a href="#">O43155</a>
Summary:	This gene encodes a member of the fibronectin leucine rich transmembrane (FLRT) family of cell adhesion molecules, which regulate early embryonic vascular and neural development. The encoded type I transmembrane protein has an extracellular region consisting of an N-terminal leucine-rich repeat domain and a type 3 fibronectin domain, followed by a transmembrane domain and a short C-terminal cytoplasmic tail domain. It functions as both a homophilic cell adhesion molecule and a heterophilic chemorepellent through its interaction with members of the uncoordinated-5 receptor family. Proteolytic removal of the extracellular region controls the migration of neurons in the developing cortex. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).