

## Product datasheet for **TL313042**

### FBXL7 Human shRNA Plasmid Kit (Locus ID 23194)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	FBXL7 Human shRNA Plasmid Kit (Locus ID 23194)
Locus ID:	23194
Synonyms:	FBL6; FBL7
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	FBXL7 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 23194). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001278317</a> , <a href="#">NM_012304</a> , <a href="#">NM_012304.1</a> , <a href="#">NM_012304.2</a> , <a href="#">NM_012304.3</a> , <a href="#">NM_012304.4</a> , <a href="#">NM_001278317.1</a> , <a href="#">BC075061</a> , <a href="#">BC075061.2</a> , <a href="#">BC050571</a> , <a href="#">BC136424</a> , <a href="#">NM_012304.5</a>
UniProt ID:	<a href="#">Q9UJT9</a>
Summary:	This gene encodes a member of the F-box protein family which is characterized by a 42-48 amino acid motif, the F-box, which binds to the S-phase kinase-associated protein 1 (Skp1) protein. The F-box proteins constitute one of the four subunits of E3 ubiquitin protein ligases called SCFs (SKP1-Cul1-F-box), which play a role in phosphorylation-dependent ubiquitination of proteins. The F-box proteins are divided into 3 subfamilies based on the other domain in the protein: F-box proteins that also have a WD-40 domain (Fbws subfamily), F-box proteins that also have leucine-rich repeats (Fbls subfamily) and F-box proteins that contain other motifs or lack known protein-interaction domains (Fbxs subfamily). The protein encoded by this gene belongs to the Fbls subfamily. Alternative splicing results in multiple transcript variants of this gene. [provided by RefSeq, May 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).