

Product datasheet for **TL313051**

Fibulin 5 (FBLN5) Human shRNA Plasmid Kit (Locus ID 10516)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Fibulin 5 (FBLN5) Human shRNA Plasmid Kit (Locus ID 10516)
Locus ID:	10516
Synonyms:	ADCL2; ARCL1A; ARMD3; DANCE; EVEC; FIBL-5; HNARMD; UP50
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	FBLN5 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 10516). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_006329 , NM_006329.1 , NM_006329.2 , NM_006329.3 , BC022280 , BC022280.1
UniProt ID:	Q9UBX5
Summary:	The protein encoded by this gene is a secreted, extracellular matrix protein containing an Arg-Gly-Asp (RGD) motif and calcium-binding EGF-like domains. It promotes adhesion of endothelial cells through interaction of integrins and the RGD motif. It is prominently expressed in developing arteries but less so in adult vessels. However, its expression is reinduced in balloon-injured vessels and atherosclerotic lesions, notably in intimal vascular smooth muscle cells and endothelial cells. Therefore, the protein encoded by this gene may play a role in vascular development and remodeling. Defects in this gene are a cause of autosomal dominant cutis laxa, autosomal recessive cutis laxa type I (CL type I), and age-related macular degeneration type 3 (ARMD3). [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL313051A virus into HEK293 cells. TL313051A virus was prepared using lenti-shRNA TL313051A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL313051B virus into HEK293 cells. TL313051B virus was prepared using lenti-shRNA TL313051B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL313051C] virus into HEK293 cells. [TL313051C] virus was prepared using lenti-shRNA [TL313051C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL313051D] virus into HEK293 cells. [TL313051D] virus was prepared using lenti-shRNA [TL313051D] and [TR30037] packaging kit.