

Product datasheet for **TL313331**

DC2L1 (DYNC2LI1) Human shRNA Plasmid Kit (Locus ID 51626)

Product data:

Product Type:	shRNA Plasmids
Product Name:	DC2L1 (DYNC2LI1) Human shRNA Plasmid Kit (Locus ID 51626)
Locus ID:	51626
Synonyms:	CGI-60; D2LIC; DKFZp564A033; dynein, cytoplasmic 2, light intermediate chain 1; dynein 2 light intermediate chain; LIC3
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	DYNC2LI1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 51626). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001012665 , NM_001193464 , NM_015522 , NM_016008 , NM_001348912 , NM_001348913 , NM_016008.1 , NM_016008.2 , NM_016008.3 , NM_015522.1 , NM_015522.2 , NM_015522.3 , NM_001012665.1 , NM_001193464.1 , BC006969 , BC006969.1 , BC016324 , BC016324.1 , BC035126 , BC040558 , BC058823 , NM_016008.4 , NM_001193464.2
UniProt ID:	Q8TCX1
Summary:	This gene encodes a protein that is a component of the dynein-2 microtubule motor protein complex that plays a role in the retrograde transport of cargo in primary cilia via the intraflagellar transport system. This gene is ubiquitously expressed and its protein, which localizes to the axoneme and Golgi apparatus, interacts directly with the cytoplasmic dynein 2 heavy chain 1 protein to form part of the multi-protein dynein-2 complex. Mutations in this gene produce defects in the dynein-2 complex which result in several types of ciliopathy including short-rib thoracic dysplasia 15 with polydactyly (SRTD15). Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Feb 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).