

## Product datasheet for **TL313342**

### DUT Human shRNA Plasmid Kit (Locus ID 1854)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	DUT Human shRNA Plasmid Kit (Locus ID 1854)
Locus ID:	1854
Synonyms:	dUTPase
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	DUT - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 1854). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001025248</a> , <a href="#">NM_001025249</a> , <a href="#">NM_001330286</a> , <a href="#">NM_001948</a> , <a href="#">NM_001025248.1</a> , <a href="#">NM_001025249.1</a> , <a href="#">NM_001948.1</a> , <a href="#">NM_001948.2</a> , <a href="#">NM_001948.3</a> , <a href="#">BC070339</a> , <a href="#">BC070339.1</a> , <a href="#">BC110377</a> , <a href="#">BC110377.1</a> , <a href="#">BC033645</a> , <a href="#">NM_001025248.2</a> , <a href="#">NM_001948.4</a>
UniProt ID:	<a href="#">P33316</a>
Summary:	This gene encodes an essential enzyme of nucleotide metabolism. The encoded protein forms a ubiquitous, homotetrameric enzyme that hydrolyzes dUTP to dUMP and pyrophosphate. This reaction serves two cellular purposes: providing a precursor (dUMP) for the synthesis of thymine nucleotides needed for DNA replication, and limiting intracellular pools of dUTP. Elevated levels of dUTP lead to increased incorporation of uracil into DNA, which induces extensive excision repair mediated by uracil glycosylase. This repair process, resulting in the removal and reincorporation of dUTP, is self-defeating and leads to DNA fragmentation and cell death. Alternative splicing of this gene leads to different isoforms that localize to either the mitochondrion or nucleus. A related pseudogene is located on chromosome 19. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).