

Product datasheet for **TL313444V**

DMC1 Human shRNA Lentiviral Particle (Locus ID 11144)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	DMC1 Human shRNA Lentiviral Particle (Locus ID 11144)
Locus ID:	11144
Synonyms:	dj199H16.1; DMC1H; LIM15
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	DMC1 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_001278208 , NM_007068 , NR_103474 , NM_007068.1 , NM_007068.2 , NM_007068.3 , NM_001278208.1 , BC125164 , BC035658 , BC125163 , BM545092 , BM843634 , NM_001363017 , NM_007068.4
UniProt ID:	Q14565
Summary:	This gene encodes a member of the superfamily of recombinases (also called DNA strand-exchange proteins). Recombinases are important for repairing double-strand DNA breaks during mitosis and meiosis. This protein, which is evolutionarily conserved, is reported to be essential for meiotic homologous recombination and may thus play an important role in generating diversity of genetic information. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).