

## Product datasheet for **TL313515**

### DDX5 Human shRNA Plasmid Kit (Locus ID 1655)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	DDX5 Human shRNA Plasmid Kit (Locus ID 1655)
Locus ID:	1655
Synonyms:	G17P1; HLR1; HUMP68; p68
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	DDX5 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 1655). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001320595</a> , <a href="#">NM_001320596</a> , <a href="#">NM_001320597</a> , <a href="#">NM_004396</a> , <a href="#">NM_004396.1</a> , <a href="#">NM_004396.2</a> , <a href="#">NM_004396.3</a> , <a href="#">NM_004396.4</a> , <a href="#">BC016027</a> , <a href="#">BC016027.1</a> , <a href="#">BM695261</a> , <a href="#">NM_004396.5</a>
UniProt ID:	<a href="#">P17844</a>
Summary:	This gene encodes a member of the DEAD box family of RNA helicases that are involved in a variety of cellular processes as a result of its role as an adaptor molecule, promoting interactions with a large number of other factors. This protein is involved in pathways that include the alteration of RNA structures, plays a role as a coregulator of transcription, a regulator of splicing, and in the processing of small noncoding RNAs. Members of this family contain nine conserved motifs, including the conserved Asp-Glu-Ala-Asp (DEAD) motif, important to ATP binding and hydrolysis as well as RNA binding and unwinding activities. Dysregulation of this gene may play a role in cancer development. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).