

## Product datasheet for **TL313639V**

### Cullin 4B (CUL4B) Human shRNA Lentiviral Particle (Locus ID 8450)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Cullin 4B (CUL4B) Human shRNA Lentiviral Particle (Locus ID 8450)
Locus ID:	8450
Synonyms:	CUL-4B; MRXHF2; MRXS15; MRXSC; SFM2
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	CUL4B - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_001079872</a> , <a href="#">NM_003588</a> , <a href="#">NM_001330624</a> , <a href="#">NM_003588.1</a> , <a href="#">NM_003588.2</a> , <a href="#">NM_003588.3</a> , <a href="#">NM_001079872.1</a> , <a href="#">BC036216</a> , <a href="#">BC036216.1</a> , <a href="#">NM_001369145</a> , <a href="#">NM_001079872.2</a>
UniProt ID:	<a href="#">Q13620</a>
Summary:	This gene is a member of the cullin family. The encoded protein forms a complex that functions as an E3 ubiquitin ligase and catalyzes the polyubiquitination of specific protein substrates in the cell. The protein interacts with a ring finger protein, and is required for the proteolysis of several regulators of DNA replication including chromatin licensing and DNA replication factor 1 and cyclin E. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).