

Product datasheet for **TL313679V**

CSTF3 Human shRNA Lentiviral Particle (Locus ID 1479)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	CSTF3 Human shRNA Lentiviral Particle (Locus ID 1479)
Locus ID:	1479
Synonyms:	CSTF-77
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	CSTF3 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	<u>NM_001033505</u> , <u>NM_001033506</u> , <u>NM_001326</u> , <u>NM_001326.1</u> , <u>NM_001326.2</u> , <u>NM_001033506.1</u> , <u>NM_001033505.1</u> , <u>BC108319</u> , <u>BC009792</u> , <u>BC010533</u> , <u>BC059948</u> , <u>NM_001326.3</u>
UniProt ID:	<u>Q12996</u>
Summary:	The protein encoded by this gene is one of three (including CSTF1 and CSTF2) cleavage stimulation factors that combine to form the cleavage stimulation factor complex (CSTF). This complex is involved in the polyadenylation and 3' end cleavage of pre-mRNAs. The encoded protein functions as a homodimer and interacts directly with both CSTF1 and CSTF2 in the CSTF complex. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).