

## Product datasheet for **TL313696V**

### Casein Kinase 2 beta (CSNK2B) Human shRNA Lentiviral Particle (Locus ID 1460)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Casein Kinase 2 beta (CSNK2B) Human shRNA Lentiviral Particle (Locus ID 1460)
Locus ID:	1460
Synonyms:	CK2B; CK2N; Ckb1; Ckb2; CSK2B; G5A; POBINDS
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	CSNK2B - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_001282385</a> , <a href="#">NM_001320</a> , <a href="#">NM_001320.1</a> , <a href="#">NM_001320.2</a> , <a href="#">NM_001320.3</a> , <a href="#">NM_001320.4</a> , <a href="#">NM_001320.6</a> , <a href="#">NM_001282385.1</a> , <a href="#">BC112017</a> , <a href="#">BC112017.1</a> , <a href="#">BC112019</a> , <a href="#">NM_001320.7</a>
UniProt ID:	<a href="#">P67870</a>
Summary:	This gene encodes the beta subunit of casein kinase II, a ubiquitous protein kinase which regulates metabolic pathways, signal transduction, transcription, translation, and replication. The enzyme is composed of three subunits, alpha, alpha prime and beta, which form a tetrameric holoenzyme. The alpha and alpha prime subunits are catalytic, while the beta subunit serves regulatory functions. The enzyme localizes to the endoplasmic reticulum and the Golgi apparatus. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .




[View online »](#)


**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

**Product images:**


GFP signal was observed under microscope at 48 hours after transduction of TL313696A virus into HEK293 cells. TL313696A virus was prepared using lenti-shRNA TL313696A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL313696B virus into HEK293 cells. TL313696B virus was prepared using lenti-shRNA TL313696B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL313696C] virus into HEK293 cells. [TL313696C] virus was prepared using lenti-shRNA [TL313696C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL313696D] virus into HEK293 cells. [TL313696D] virus was prepared using lenti-shRNA [TL313696D] and [TR30037] packaging kit.