

Product datasheet for **TL313741**

CRAT Human shRNA Plasmid Kit (Locus ID 1384)

Product data:

Product Type:	shRNA Plasmids
Product Name:	CRAT Human shRNA Plasmid Kit (Locus ID 1384)
Locus ID:	1384
Synonyms:	CAT; CAT1; NBIA8
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	CRAT - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 1384). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_000755 , NM_001257363 , NM_004003 , NM_144782 , NR_028048 , NM_001346546 , NM_001346547 , NM_001346548 , NM_001346549 , NM_000755.1 , NM_000755.2 , NM_000755.3 , NM_000755.4 , NM_004003.1 , NM_004003.2 , NM_004003.3 , NM_001257363.1 , NM_001257363.2 , NM_144782.1 , BC000723 , BM553197 , NM_000755.5
UniProt ID:	P43155
Summary:	This gene encodes carnitine O-acetyltransferase, a member of the carnitine acyltransferase family and a key metabolic pathway enzyme which plays an important role in energy homeostasis and fat metabolism. This enzyme catalyzes the reversible transfer of acyl groups from an acyl-CoA thioester to carnitine and regulates the ratio of acyl-CoA/CoA. It is found in both the mitochondria and the peroxisome. Alternative splicing results in transcript variants encoding different isoforms that may localize to different subcellular compartments. [provided by RefSeq, Oct 2016]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).