

Product datasheet for **TL313927**

CHRM2 Human shRNA Plasmid Kit (Locus ID 1129)

Product data:

Product Type:	shRNA Plasmids
Product Name:	CHRM2 Human shRNA Plasmid Kit (Locus ID 1129)
Locus ID:	1129
Synonyms:	HM2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	CHRM2 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 1129). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<u>NM_000739</u> , <u>NM_001006626</u> , <u>NM_001006627</u> , <u>NM_001006628</u> , <u>NM_001006629</u> , <u>NM_001006630</u> , <u>NM_001006631</u> , <u>NM_001006632</u> , <u>NM_001006633</u> , <u>NM_001006629.1</u> , <u>NM_000739.1</u> , <u>NM_000739.2</u> , <u>NM_001006628.1</u> , <u>NM_001006631.1</u> , <u>NM_001006627.1</u> , <u>NM_001006630.1</u> , <u>NM_001006626.1</u> , <u>NM_001006633.1</u> , <u>BC095547</u> , <u>BC106741</u> , <u>BC106742</u> , <u>NM_001006629.3</u> , <u>NM_001006632.3</u> , <u>NM_001006631.3</u> , <u>NM_001006628.3</u> , <u>NM_001006627.3</u> , <u>NM_000739.3</u> , <u>NM_001006626.3</u>
UniProt ID:	<u>P08172</u>
Summary:	The muscarinic cholinergic receptors belong to a larger family of G protein-coupled receptors. The functional diversity of these receptors is defined by the binding of acetylcholine to these receptors and includes cellular responses such as adenylylase inhibition, phosphoinositide degeneration, and potassium channel mediation. Muscarinic receptors influence many effects of acetylcholine in the central and peripheral nervous system. The muscarinic cholinergic receptor 2 is involved in mediation of bradycardia and a decrease in cardiac contractility. Multiple alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Jul 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).