

## Product datasheet for **TL314122**

### CCR9 Human shRNA Plasmid Kit (Locus ID 10803)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	CCR9 Human shRNA Plasmid Kit (Locus ID 10803)
Locus ID:	10803
Synonyms:	CC-CKR-9; CDw199; GPR-9-6; GPR28
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	CCR9 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 10803). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001256369</a> , <a href="#">NM_006641</a> , <a href="#">NM_031200</a> , <a href="#">NR_036525</a> , <a href="#">NM_006641.1</a> , <a href="#">NM_006641.2</a> , <a href="#">NM_006641.3</a> , <a href="#">NM_031200.1</a> , <a href="#">NM_031200.2</a> , <a href="#">NM_001256369.1</a> , <a href="#">BC069678</a> , <a href="#">BC095516</a> , <a href="#">NM_031200.3</a>
UniProt ID:	<a href="#">P51686</a>
Summary:	The protein encoded by this gene is a G protein-coupled receptor with seven transmembrane domains that belongs to the beta chemokine receptor family. Chemokines and their receptors are key regulators of thymocyte migration and maturation in normal and inflammation conditions. This gene is differentially expressed in T lymphocytes of the small intestine and colon, and its interaction with chemokine 25 contributes to intestinal intra-epithelial lymphocyte homing to the small intestine. This suggests a role for this gene in directing immune responses to different segments of the gastrointestinal tract. This gene and its exclusive ligand, chemokine 25, are overexpressed in a variety of malignant tumors and are closely associated with tumor proliferation, apoptosis, invasion, migration and drug resistance. This gene maps to the chemokine receptor gene cluster. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2020]



[View online »](#)

- shRNA Design:** These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact [techsupport@origene.com](mailto:techsupport@origene.com). If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).
- Performance Guaranteed:** OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.
- For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).