

Product datasheet for **TL314194**

CARD8 Human shRNA Plasmid Kit (Locus ID 22900)

Product data:

Product Type:	shRNA Plasmids
Product Name:	CARD8 Human shRNA Plasmid Kit (Locus ID 22900)
Locus ID:	22900
Synonyms:	CARDINAL; DACAR; DAKAR; NDPP; NDPP1; TUCAN
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	CARD8 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 22900). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001184900 , NM_001184901 , NM_001184902 , NM_001184903 , NM_001184904 , NM_014959 , NR_033678 , NR_033679 , NR_033680 , NM_001351782 , NM_001351783 , NM_001351784 , NM_001351786 , NM_001351787 , NM_001351788 , NM_001351789 , NM_001351790 , NM_001351791 , NM_001351792 , NR_147707 , NM_014959.1 , NM_014959.2 , NM_014959.3 , NM_001184904.1 , NM_001184902.1 , NM_001184903.1 , NM_001184901.1 , NM_001184900.1 , BC056891 , BC045539 , BC034564 , BC041697 , BC046136 , NM_001365950 , NM_001184900.3 , NM_014959.5 , NM_001184902.2
UniProt ID:	Q9Y2G2
Summary:	The protein encoded by this gene belongs to the caspase recruitment domain (CARD)-containing family of proteins, which are involved in pathways leading to activation of caspases or nuclear factor kappa-B (NFkB). This protein may be a component of the inflammasome, a protein complex that plays a role in the activation of proinflammatory caspases. It is thought that this protein acts as an adaptor molecule that negatively regulates NFkB activation, CASP1-dependent IL1B secretion, and apoptosis. Polymorphisms in this gene may be associated with a susceptibility to rheumatoid arthritis. Alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, May 2010]



[View online »](#)

shRNA Design:

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).