

Product datasheet for TL314543

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

beta 2 Microglobulin (B2M) Human shRNA Plasmid Kit (Locus ID 567)

Product data:

Product Type: shRNA Plasmids

Product Name: beta 2 Microglobulin (B2M) Human shRNA Plasmid Kit (Locus ID 567)

Locus ID: 567

Synonyms: IMD43

Vector: pGFP-C-shLenti (TR30023)

E. coli Selection: Chloramphenicol (34 ug/ml)

Mammalian Cell

Selection:

Puromycin

Format: Lentiviral plasmids

Components: B2M - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 567). 5µg

purified plasmid DNA per construct

29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.

RefSeq: NM 004048, NM 004048.1, NM 004048.2, BC032589, BC032589.1, BC064910, NM 004048.4

UniProt ID: P61769

Summary: This gene encodes a serum protein found in association with the major histocompatibility

complex (MHC) class I heavy chain on the surface of nearly all nucleated cells. The protein has

a predominantly beta-pleated sheet structure that can form amyloid fibrils in some

pathological conditions. The encoded antimicrobial protein displays antibacterial activity in

amniotic fluid. A mutation in this gene has been shown to result in hypercatabolic

hypoproteinemia.[provided by RefSeq, Aug 2014]

shRNA Design: These shRNA constructs were designed against multiple splice variants at this gene locus. To

be certain that your variant of interest is targeted, please contact <u>techsupport@origene.com</u>. If you need a special design or shRNA sequence, please utilize our <u>custom shRNA service</u>.






Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:




GFP signal was observed under microscope at 48 hours after transduction of TL314543A virus into HEK293 cells. TL314543A virus was prepared using lenti-shRNA TL314543A and [TR30037] packaging kit.






GFP signal was observed under microscope at 48 hours after transduction of TL314543B virus into HEK293 cells. TL314543B virus was prepared using lenti-shRNA TL314543B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL314543C] virus into HEK293 cells. [TL314543C] virus was prepared using lenti-shRNA [TL314543C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL314543D] virus into HEK293 cells. [TL314543D] virus was prepared using lenti-shRNA [TL314543D] and [TR30037] packaging kit.