

## Product datasheet for **TL314596**

### ATP1B4 Human shRNA Plasmid Kit (Locus ID 23439)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	ATP1B4 Human shRNA Plasmid Kit (Locus ID 23439)
Locus ID:	23439
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	ATP1B4 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 23439). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001142447</a> , <a href="#">NM_012069</a> , <a href="#">NM_012069.1</a> , <a href="#">NM_012069.2</a> , <a href="#">NM_012069.3</a> , <a href="#">NM_012069.4</a> , <a href="#">NM_001142447.1</a> , <a href="#">NM_001142447.2</a> , <a href="#">BC117227</a> , <a href="#">BC143404</a> , <a href="#">BC143406</a> , <a href="#">BC143407</a> , <a href="#">NM_001142447.3</a> , <a href="#">NM_012069.5</a>
UniProt ID:	<a href="#">Q9UN42</a>
Summary:	This gene has been found in all vertebrate genomes sequenced to date. However, this gene has undergone a change in function in placental mammals compared to other species. Specifically, in fish, avian, and amphibian species, this gene encodes plasma membrane-bound beta-subunits of Na,K-ATPase. In placental mammals, the encoded protein interacts with the nuclear transcriptional coregulator SKIP and may be involved in the regulation of TGF-beta signaling. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).