

Product datasheet for **TL314703**

Amphiregulin (AREG) Human shRNA Plasmid Kit (Locus ID 374)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Amphiregulin (AREG) Human shRNA Plasmid Kit (Locus ID 374)
Locus ID:	374
Synonyms:	AR; AREGB; CRDGF; SDGF
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	AREG - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 374). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001657 , NM_001657.1 , NM_001657.2 , NM_001657.3 , BC009799 , BC146953 , BC146967 , NM_001657.4
UniProt ID:	P15514
Summary:	The protein encoded by this gene is a member of the epidermal growth factor family. It is an autocrine growth factor as well as a mitogen for astrocytes, Schwann cells and fibroblasts. It is related to epidermal growth factor (EGF) and transforming growth factor alpha (TGF-alpha). The protein interacts with the EGF/TGF-alpha receptor to promote the growth of normal epithelial cells, and it inhibits the growth of certain aggressive carcinoma cell lines. It also functions in mammary gland, oocyte and bone tissue development. This gene is associated with a psoriasis-like skin phenotype, and is also associated with other pathological disorders, including various types of cancers and inflammatory conditions. [provided by RefSeq, Apr 2014]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL314703A virus into HEK293 cells. TL314703A virus was prepared using lenti-shRNA TL314703A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL314703B virus into HEK293 cells. TL314703B virus was prepared using lenti-shRNA TL314703B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL314703C] virus into HEK293 cells. [TL314703C] virus was prepared using lenti-shRNA [TL314703C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL314703D] virus into HEK293 cells. [TL314703D] virus was prepared using lenti-shRNA [TL314703D] and [TR30037] packaging kit.