

Product datasheet for **TL314731**

APLP2 Human shRNA Plasmid Kit (Locus ID 334)

Product data:

Product Type:	shRNA Plasmids
Product Name:	APLP2 Human shRNA Plasmid Kit (Locus ID 334)
Locus ID:	334
Synonyms:	APLP-2; APPH; APPL2; CDEBP
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	APLP2 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 334). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC004371 , NM_001142276 , NM_001142277 , NM_001142278 , NM_001243299 , NM_001642 , NR_024515 , NR_024516 , NM_001328682 , NM_001328684 , NM_001328685 , NM_001328686 , NM_001642.1 , NM_001642.2 , NM_001142276.1 , NM_001142277.1 , NM_001142278.1 , NM_001243299.1 , BC004371.1 , BC000373 , BM992563 , NM_001142278.2 , NM_001642.3 , NM_001142277.2 , NM_001142276.2
UniProt ID:	Q06481
Summary:	This gene encodes amyloid precursor- like protein 2 (APLP2), which is a member of the APP (amyloid precursor protein) family including APP, APLP1 and APLP2. This protein is ubiquitously expressed. It contains heparin-, copper- and zinc- binding domains at the N-terminus, BPTI/Kunitz inhibitor and E2 domains in the middle region, and transmembrane and intracellular domains at the C-terminus. This protein interacts with major histocompatibility complex (MHC) class I molecules. The synergy of this protein and the APP is required to mediate neuromuscular transmission, spatial learning and synaptic plasticity. This protein has been implicated in the pathogenesis of Alzheimer's disease. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Aug 2011]



[View online »](#)

shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .
Performance Guaranteed:	<p>OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.</p> <p>For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).</p>