

Product datasheet for **TL314932**

ADFP (PLIN2) Human shRNA Plasmid Kit (Locus ID 123)

Product data:

Product Type:	shRNA Plasmids
Product Name:	ADFP (PLIN2) Human shRNA Plasmid Kit (Locus ID 123)
Locus ID:	123
Synonyms:	ADFP; ADRP
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	PLIN2 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 123). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001122 , NR_038064 , NM_001122.1 , NM_001122.2 , NM_001122.3 , BC005127 , BC005127.2 , NM_001122.4
UniProt ID:	Q99541
Summary:	The protein encoded by this gene belongs to the perilipin family, members of which coat intracellular lipid storage droplets. This protein is associated with the lipid globule surface membrane material, and maybe involved in development and maintenance of adipose tissue. However, it is not restricted to adipocytes as previously thought, but is found in a wide range of cultured cell lines, including fibroblasts, endothelial and epithelial cells, and tissues, such as lactating mammary gland, adrenal cortex, Sertoli and Leydig cells, and hepatocytes in alcoholic liver cirrhosis, suggesting that it may serve as a marker of lipid accumulation in diverse cell types and diseases. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Mar 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL314932A virus into HEK293 cells. TL314932A virus was prepared using lenti-shRNA TL314932A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL314932B virus into HEK293 cells. TL314932B virus was prepared using lenti-shRNA TL314932B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL314932C] virus into HEK293 cells. [TL314932C] virus was prepared using lenti-shRNA [TL314932C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL314932D] virus into HEK293 cells. [TL314932D] virus was prepared using lenti-shRNA [TL314932D] and [TR30037] packaging kit.