

Product datasheet for **TL314976V**

beta Actin (ACTB) Human shRNA Lentiviral Particle (Locus ID 60)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	beta Actin (ACTB) Human shRNA Lentiviral Particle (Locus ID 60)
Locus ID:	60
Synonyms:	BRWS1; PS1TP5BP1
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	ACTB - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_001101 , NM_001101.1 , NM_001101.2 , NM_001101.3 , BC013380 , BC013380.2 , BC001301 , BC002409 , BC004251 , BC008633 , BC009636 , BC012854 , BC013835 , BC014401 , BC014861 , BC016045 , BC023204 , BC113036 , NM_001101.5
UniProt ID:	P60709
Summary:	This gene encodes one of six different actin proteins. Actins are highly conserved proteins that are involved in cell motility, structure, integrity, and intercellular signaling. The encoded protein is a major constituent of the contractile apparatus and one of the two nonmuscle cytoskeletal actins that are ubiquitously expressed. Mutations in this gene cause Baraitser-Winter syndrome 1, which is characterized by intellectual disability with a distinctive facial appearance in human patients. Numerous pseudogenes of this gene have been identified throughout the human genome. [provided by RefSeq, Aug 2017]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).