

Product datasheet for **TL314983**

Prostatic Acid Phosphatase (ACPP) Human shRNA Plasmid Kit (Locus ID 55)

Product data:

Product Type:	shRNA Plasmids
Product Name:	Prostatic Acid Phosphatase (ACPP) Human shRNA Plasmid Kit (Locus ID 55)
Locus ID:	55
Synonyms:	5'-NT; ACP-3; ACP3
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	ACPP - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 55). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001099 , NM_001134194 , NM_001292037 , NM_001099.1 , NM_001099.2 , NM_001099.3 , NM_001099.4 , NM_001134194.1 , NM_001292037.1 , BC016344 , BC016344.1 , BC007460 , BC008493 , BC017877 , NM_001292037.2 , NM_001134194.2 , NM_001099.5
UniProt ID:	P15309
Summary:	This gene encodes an enzyme that catalyzes the conversion of orthophosphoric monoester to alcohol and orthophosphate. It is synthesized under androgen regulation and is secreted by the epithelial cells of the prostate gland. An alternatively spliced transcript variant encoding a longer isoform has been found for this gene. This isoform contains a transmembrane domain and is localized in the plasma membrane-endosomal-lysosomal pathway. [provided by RefSeq, Sep 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).