

Product datasheet for **TL315028V**

ABCB7 Human shRNA Lentiviral Particle (Locus ID 22)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	ABCB7 Human shRNA Lentiviral Particle (Locus ID 22)
Locus ID:	22
Synonyms:	ABC7; ASAT; Atm1p; EST140535
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	ABCB7 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	<u>NM_001271696</u> , <u>NM_001271697</u> , <u>NM_001271698</u> , <u>NM_001271699</u> , <u>NM_004299</u> , <u>NM_004299.1</u> , <u>NM_004299.2</u> , <u>NM_004299.3</u> , <u>NM_004299.4</u> , <u>NM_001271697.1</u> , <u>NM_001271698.1</u> , <u>NM_001271696.1</u> , <u>NM_001271699.1</u> , <u>BC006323</u> , <u>BC006323.2</u> , <u>NM_001271697.3</u> , <u>NM_001271698.3</u> , <u>NM_001271699.3</u> , <u>NM_004299.6</u> , <u>NM_001271696.3</u>
UniProt ID:	<u>O75027</u>
Summary:	The membrane-associated protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the MDR/TAP subfamily. Members of the MDR/TAP subfamily are involved in multidrug resistance as well as antigen presentation. This gene encodes a half-transporter involved in the transport of heme from the mitochondria to the cytosol. With iron/sulfur cluster precursors as its substrates, this protein may play a role in metal homeostasis. Mutations in this gene have been associated with mitochondrial iron accumulation and isodicentric (X)(q13) and sideroblastic anemia. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Nov 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).