

## Product datasheet for **TL315126**

### EPHA10 Human shRNA Plasmid Kit (Locus ID 284656)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	EPHA10 Human shRNA Plasmid Kit (Locus ID 284656)
Locus ID:	284656
Synonyms:	EphA10s protein; EPH receptor A10; Ephrin type-A receptor 10; FLJ16103; FLJ16103, FLJ33655, MGC43817; FLJ33655; MGC43817; OTTHUMP00000004519
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	EPHA10 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 284656). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001004338</a> , <a href="#">NM_001099439</a> , <a href="#">NM_173641</a> , <a href="#">NM_173641.1</a> , <a href="#">NM_173641.2</a> , <a href="#">NM_001099439.1</a> , <a href="#">NM_001004338.1</a> , <a href="#">BC067734</a> , <a href="#">BC067734.1</a> , <a href="#">BC039363</a> , <a href="#">BC041360</a> , <a href="#">BC112933</a> , <a href="#">NM_001099439.2</a>
UniProt ID:	<a href="#">Q5JZY3</a>
Summary:	Ephrin receptors, the largest subfamily of receptor tyrosine kinases (RTKs), and their ephrin ligands are important mediators of cell-cell communication regulating cell attachment, shape, and mobility in neuronal and epithelial cells (Aasheim et al., 2005 [PubMed 15777695]). See MIM 179610 for additional background on Eph receptors and ephrins.[supplied by OMIM, Mar 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).