

Product datasheet for **TL315501**

MUM1 (IRF4) Human shRNA Plasmid Kit (Locus ID 3662)

Product data:

Product Type:	shRNA Plasmids
Product Name:	MUM1 (IRF4) Human shRNA Plasmid Kit (Locus ID 3662)
Locus ID:	3662
Synonyms:	LSIRF; MUM1; NF-EM5; SHEP8
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	IRF4 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 3662). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	BC015752 , NM_001195286 , NM_002460 , NR_036585 , NR_046000 , NM_002460.1 , NM_002460.2 , NM_002460.3 , NM_001195286.1 , BC015752.1 , BM551611 , NM_001195286.2 , NM_002460.4
UniProt ID:	Q15306
Summary:	The protein encoded by this gene belongs to the IRF (interferon regulatory factor) family of transcription factors, characterized by an unique tryptophan pentad repeat DNA-binding domain. The IRFs are important in the regulation of interferons in response to infection by virus, and in the regulation of interferon-inducible genes. This family member is lymphocyte specific and negatively regulates Toll-like-receptor (TLR) signaling that is central to the activation of innate and adaptive immune systems. A chromosomal translocation involving this gene and the IgH locus, t(6;14)(p25;q32), may be a cause of multiple myeloma. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Aug 2010]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .




[View online »](#)


**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

Product images:


GFP signal was observed under microscope at 48 hours after transduction of TL315501A virus into HEK293 cells. TL315501A virus was prepared using lenti-shRNA TL315501A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL315501B virus into HEK293 cells. TL315501B virus was prepared using lenti-shRNA TL315501B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL315501C] virus into HEK293 cells. [TL315501C] virus was prepared using lenti-shRNA [TL315501C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL315501D] virus into HEK293 cells. [TL315501D] virus was prepared using lenti-shRNA [TL315501D] and [TR30037] packaging kit.