

Product datasheet for **TL315609**

SRCAP Human shRNA Plasmid Kit (Locus ID 10847)

Product data:

Product Type:	shRNA Plasmids
Product Name:	SRCAP Human shRNA Plasmid Kit (Locus ID 10847)
Locus ID:	10847
Synonyms:	DOMO1; EAF1; FLHS; SWR1
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	SRCAP - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 10847). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_006662 , NM_006662.1 , NM_006662.2 , BC025284 , BC063578 , BC172428 , BM677852
UniProt ID:	Q6ZRS2
Summary:	This gene encodes the core catalytic component of the multiprotein chromatin-remodeling SRCAP complex. The encoded protein is an ATPase that is necessary for the incorporation of the histone variant H2A.Z into nucleosomes. It can function as a transcriptional activator in Notch-mediated, CREB-mediated and steroid receptor-mediated transcription. Mutations in this gene cause Floating-Harbor syndrome, a rare disorder characterized by short stature, language deficits and dysmorphic facial features. [provided by RefSeq, Feb 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).