

## Product datasheet for **TL315696**

### **KMT2B Human shRNA Plasmid Kit (Locus ID 9757)**

#### **Product data:**

<b>Product Type:</b>	shRNA Plasmids
<b>Product Name:</b>	KMT2B Human shRNA Plasmid Kit (Locus ID 9757)
<b>Locus ID:</b>	9757
<b>Synonyms:</b>	CXXC10; DYT28; HRX2; MLL1B; MLL2; MLL4; TRX2; WBP-7; WBP7
<b>Vector:</b>	pGFP-C-shLenti (TR30023)
<b>E. coli Selection:</b>	Chloramphenicol (34 ug/ml)
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Format:</b>	Lentiviral plasmids
<b>Components:</b>	KMT2B - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 9757). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
<b>RefSeq:</b>	<a href="#">NM_014727</a> , <a href="#">NM_014727.1</a> , <a href="#">NM_014727.2</a> , <a href="#">BC007353</a> , <a href="#">BC009337</a> , <a href="#">BC026861</a> , <a href="#">BM967067</a> , <a href="#">NM_014727.3</a>
<b>UniProt ID:</b>	<a href="#">Q9UMN6</a>
<b>Summary:</b>	This gene encodes a protein which contains multiple domains including a CXXC zinc finger, three PHD zinc fingers, two FY-rich domains, and a SET (suppressor of variegation, enhancer of zeste, and trithorax) domain. The SET domain is a conserved C-terminal domain that characterizes proteins of the MLL (mixed-lineage leukemia) family. This gene is ubiquitously expressed in adult tissues. It is also amplified in solid tumor cell lines, and may be involved in human cancer. Two alternatively spliced transcript variants encoding distinct isoforms have been reported for this gene, however, the full length nature of the shorter transcript is not known. [provided by RefSeq, Jul 2008]
<b>shRNA Design:</b>	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .




[View online »](#)


**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.


For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).

**Product images:**


GFP signal was observed under microscope at 48 hours after transduction of TL315696A virus into HEK293 cells. TL315696A virus was prepared using lenti-shRNA TL315696A and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of TL315696B virus into HEK293 cells. TL315696B virus was prepared using lenti-shRNA TL315696B and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL315696C] virus into HEK293 cells. [TL315696C] virus was prepared using lenti-shRNA [TL315696C] and [TR30037] packaging kit.



GFP signal was observed under microscope at 48 hours after transduction of [TL315696D] virus into HEK293 cells. [TL315696D] virus was prepared using lenti-shRNA [TL315696D] and [TR30037] packaging kit.