

Product datasheet for **TL316441V**

AMFR Human shRNA Lentiviral Particle (Locus ID 267)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	AMFR Human shRNA Lentiviral Particle (Locus ID 267)
Locus ID:	267
Synonyms:	GP78; RNF45
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	AMFR - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_001144 , NM_001323511 , NM_001323512 , NM_138958 , NM_001144.1 , NM_001144.2 , NM_001144.3 , NM_001144.4 , NM_001144.5 , BC069197 , BC069197.1 , BC017043 , BC051032 , BC056869 , NM_001144.6
UniProt ID:	Q9UKV5
Summary:	This locus encodes a glycosylated transmembrane receptor. Its ligand, autocrine motility factor, is a tumor motility-stimulating protein secreted by tumor cells. The encoded receptor is also a member of the E3 ubiquitin ligase family of proteins. It catalyzes ubiquitination and endoplasmic reticulum-associated degradation of specific proteins. [provided by RefSeq, Feb 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).