

## Product datasheet for **TL316453V**

### Asialoglycoprotein Receptor 1 (ASGR1) Human shRNA Lentiviral Particle (Locus ID 432)

#### Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	Asialoglycoprotein Receptor 1 (ASGR1) Human shRNA Lentiviral Particle (Locus ID 432)
Locus ID:	432
Synonyms:	ASGPR; ASGPR1; CLEC4H1; HL-1
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	ASGR1 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 <sup>7</sup> TU/ml.
RefSeq:	<a href="#">NM_001197216</a> , <a href="#">NM_001671</a> , <a href="#">NM_001671.1</a> , <a href="#">NM_001671.2</a> , <a href="#">NM_001671.3</a> , <a href="#">NM_001671.4</a> , <a href="#">NM_001197216.1</a> , <a href="#">NM_001197216.2</a> , <a href="#">BC032130</a> , <a href="#">BC032130.2</a> , <a href="#">BM818592</a> , <a href="#">NM_001671.5</a>
UniProt ID:	<a href="#">P07306</a>
Summary:	This gene encodes a subunit of the asialoglycoprotein receptor. This receptor is a transmembrane protein that plays a critical role in serum glycoprotein homeostasis by mediating the endocytosis and lysosomal degradation of glycoproteins with exposed terminal galactose or N-acetylgalactosamine residues. The asialoglycoprotein receptor may facilitate hepatic infection by multiple viruses including hepatitis B, and is also a target for liver-specific drug delivery. The asialoglycoprotein receptor is a hetero-oligomeric protein composed of major and minor subunits, which are encoded by different genes. The protein encoded by this gene is the more abundant major subunit. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jan 2011]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).