

## Product datasheet for **TL316509**

### Cellular Apoptosis Susceptibility (CSE1L) Human shRNA Plasmid Kit (Locus ID 1434)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	Cellular Apoptosis Susceptibility (CSE1L) Human shRNA Plasmid Kit (Locus ID 1434)
Locus ID:	1434
Synonyms:	CAS; CSE1; XPO2
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	CSE1L - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 1434). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001256135</a> , <a href="#">NM_001316</a> , <a href="#">NM_177436</a> , <a href="#">NR_045796</a> , <a href="#">NM_001316.1</a> , <a href="#">NM_001316.2</a> , <a href="#">NM_001316.3</a> , <a href="#">NM_001256135.1</a> , <a href="#">NM_177436.1</a> , <a href="#">BC109313</a> , <a href="#">BC108309</a> , <a href="#">BC109314</a> , <a href="#">NM_001362762</a> , <a href="#">NM_001316.4</a>
UniProt ID:	<a href="#">P55060</a>
Summary:	Proteins that carry a nuclear localization signal (NLS) are transported into the nucleus by the importin-alpha/beta heterodimer. Importin-alpha binds the NLS, while importin-beta mediates translocation through the nuclear pore complex. After translocation, RanGTP binds importin-beta and displaces importin-alpha. Importin-alpha must then be returned to the cytoplasm, leaving the NLS protein behind. The protein encoded by this gene binds strongly to NLS-free importin-alpha, and this binding is released in the cytoplasm by the combined action of RANBP1 and RANGAP1. In addition, the encoded protein may play a role both in apoptosis and in cell proliferation. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jan 2012]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).