

Product datasheet for **TL316537V**

FGF1 Human shRNA Lentiviral Particle (Locus ID 2246)

Product data:

Product Type:	shRNA Lentiviral Particles
Product Name:	FGF1 Human shRNA Lentiviral Particle (Locus ID 2246)
Locus ID:	2246
Synonyms:	AFGF; ECGF; ECGF-beta; ECGFA; ECGFB; FGF-1; FGF-alpha; FGFA; GLIO703; HBGF-1; HBGF1
Vector:	pGFP-C-shLenti (TR30023)
Format:	Lentiviral particles
Components:	FGF1 - Human shRNA lentiviral particles (4 unique 29mer target-specific shRNA, 1 scramble control), 0.5 ml each, >10 ⁷ TU/ml.
RefSeq:	NM_000800 , NM_001144892 , NM_001144934 , NM_001144935 , NM_001257205 , NM_001257206 , NM_001257207 , NM_001257208 , NM_001257209 , NM_001257210 , NM_001257211 , NM_001257212 , NM_033136 , NM_033137 , NR_026695 , NR_026696 , NM_001354951 , NM_001354952 , NM_001354953 , NM_001354954 , NM_001354955 , NM_001354956 , NM_001354957 , NM_001354958 , NM_001354959 , NM_001354961 , NM_001354962 , NM_001354963 , NM_001354964 , NM_000800.1 , NM_000800.2 , NM_000800.3 , NM_000800.4 , NM_033137.1 , NM_033137.2 , NM_033136.1 , NM_033136.2 , NM_033136.3 , NM_001144892.1 , NM_001144892.2 , NM_001144934.1 , NM_001144935.1 , NM_001257206.1 , NM_001257211.1 , NM_001257212.1 , NM_001257205.1 , NM_001257207.1 , NM_001257208.1 , NM_001257209.1 , NM_001257210.1 , BC032697 , BC032697.1 , NM_001144934.2 , NM_000800.5 , NM_033136.4 , NM_001144935.2 , NM_001257208.2 , NM_001257212.2 , NM_001144892.3 , NM_033137.4 , NM_001257211.2 , NM_001257206.2
UniProt ID:	P05230
Summary:	The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. This protein functions as a modifier of endothelial cell migration and proliferation, as well as an angiogenic factor. It acts as a mitogen for a variety of mesoderm- and neuroectoderm-derived cells in vitro, thus is thought to be involved in organogenesis. Multiple alternatively spliced variants encoding different isoforms have been described. [provided by RefSeq, Jan 2009]



[View online »](#)

shRNA Design:

These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com. If you need a special design or shRNA sequence, please utilize our [custom shRNA service](#).

Performance Guaranteed:

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).