

## Product datasheet for **TL316552**

### **GATA4 Human shRNA Plasmid Kit (Locus ID 2626)**

#### **Product data:**

<b>Product Type:</b>	shRNA Plasmids
<b>Product Name:</b>	GATA4 Human shRNA Plasmid Kit (Locus ID 2626)
<b>Locus ID:</b>	2626
<b>Synonyms:</b>	ASD2; TACHD; TOF; VSD1
<b>Vector:</b>	pGFP-C-shLenti (TR30023)
<b>E. coli Selection:</b>	Chloramphenicol (34 ug/ml)
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Format:</b>	Lentiviral plasmids
<b>Components:</b>	GATA4 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 2626). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
<b>RefSeq:</b>	<a href="#">NM_001308093</a> , <a href="#">NM_001308094</a> , <a href="#">NM_002052</a> , <a href="#">NM_002052.1</a> , <a href="#">NM_002052.2</a> , <a href="#">NM_002052.3</a> , <a href="#">NM_002052.4</a> , <a href="#">BC101580</a> , <a href="#">BC101580.1</a> , <a href="#">BC033672</a> , <a href="#">BC068079</a> , <a href="#">BC105108</a> , <a href="#">BC143434</a> , <a href="#">BC143479</a> , <a href="#">BM560562</a> , <a href="#">BM768332</a> , <a href="#">NM_002052.5</a>
<b>UniProt ID:</b>	<a href="#">P43694</a>
<b>Summary:</b>	This gene encodes a member of the GATA family of zinc-finger transcription factors. Members of this family recognize the GATA motif which is present in the promoters of many genes. This protein is thought to regulate genes involved in embryogenesis and in myocardial differentiation and function, and is necessary for normal testicular development. Mutations in this gene have been associated with cardiac septal defects. Additionally, alterations in gene expression have been associated with several cancer types. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2015]
<b>shRNA Design:</b>	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).