

Product datasheet for **TL316630**

NBL1 Human shRNA Plasmid Kit (Locus ID 4681)

Product data:

Product Type:	shRNA Plasmids
Product Name:	NBL1 Human shRNA Plasmid Kit (Locus ID 4681)
Locus ID:	4681
Synonyms:	D1S1733E; DAN; DAND1; NB; NO3
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	NBL1 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 4681). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001204084 , NM_001204085 , NM_001204086 , NM_001278164 , NM_001278165 , NM_001278166 , NM_005380 , NM_182744 , NM_005380.1 , NM_005380.3 , NM_005380.4 , NM_005380.5 , NM_005380.6 , NM_005380.7 , NM_182744.1 , NM_182744.2 , NM_182744.3 , NM_001204084.1 , NM_001204084.2 , NM_001204085.1 , NM_001278164.1 , NM_001278165.1 , NM_001278166.1 , NM_001204086.1 , NM_001204086.2 , BC012037 , BC012037.1 , BM988954
UniProt ID:	P41271
Summary:	This gene product is the founding member of the evolutionarily conserved CAN (Cerberus and DAN) family of proteins, which contain a domain resembling the CTCK (C-terminal cystine knot-like) motif found in a number of signaling molecules. These proteins are secreted, and act as BMP (bone morphogenetic protein) antagonists by binding to BMPs and preventing them from interacting with their receptors. They may thus play an important role during growth and development. Alternatively spliced transcript variants have been identified for this gene. Read-through transcripts between this locus and the upstream mitochondrial inner membrane organizing system 1 gene (GeneID 440574) have been observed. [provided by RefSeq, May 2013]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .



[View online >](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).