

## Product datasheet for **TL316852**

### SLC16A3 Human shRNA Plasmid Kit (Locus ID 9123)

#### Product data:

Product Type:	shRNA Plasmids
Product Name:	SLC16A3 Human shRNA Plasmid Kit (Locus ID 9123)
Locus ID:	9123
Synonyms:	MCT-3; MCT-4; MCT 3; MCT3; MCT 4; MCT4
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	SLC16A3 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 9123). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	<a href="#">NM_001042422</a> , <a href="#">NM_001042423</a> , <a href="#">NM_001206950</a> , <a href="#">NM_001206951</a> , <a href="#">NM_001206952</a> , <a href="#">NM_004207</a> , <a href="#">NM_004207.1</a> , <a href="#">NM_004207.2</a> , <a href="#">NM_004207.3</a> , <a href="#">NM_001042422.1</a> , <a href="#">NM_001042422.2</a> , <a href="#">NM_001042423.1</a> , <a href="#">NM_001042423.2</a> , <a href="#">NM_001206950.1</a> , <a href="#">NM_001206951.1</a> , <a href="#">NM_001206952.1</a> , <a href="#">BC112267</a> , <a href="#">BC112269</a> , <a href="#">BM926587</a> , <a href="#">BM980920</a>
UniProt ID:	<a href="#">O15427</a>
Summary:	Lactic acid and pyruvate transport across plasma membranes is catalyzed by members of the proton-linked monocarboxylate transporter (MCT) family, which has been designated solute carrier family-16. Each MCT appears to have slightly different substrate and inhibitor specificities and transport kinetics, which are related to the metabolic requirements of the tissues in which it is found. The MCTs, which include MCT1 (SLC16A1; MIM 600682) and MCT2 (SLC16A7; MIM 603654), are characterized by 12 predicted transmembrane domains (Price et al., 1998 [PubMed 9425115]).[supplied by OMIM, Mar 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).