

## Product datasheet for **TL317168**

### Angiopoietin like 4 (ANGPTL4) Human shRNA Plasmid Kit (Locus ID 51129)

#### Product data:

|                           |   |
|---------------------------|---|
| Product Type:             | shRNA Plasmids  |
| Product Name:             | Angiopoietin like 4 (ANGPTL4) Human shRNA Plasmid Kit (Locus ID 51129)  |
| Locus ID:                 | 51129   |
| Synonyms:                 | ARP4; FIAF; HARP; HFARP; NL2; PGAR; pp1158; TGQTL; UNQ171   |
| Vector:                   | pGFP-C-shLenti (TR30023)  |
| E. coli Selection:        | Chloramphenicol (34 ug/ml)  |
| Mammalian Cell Selection: | Puromycin   |
| Format:                   | Lentiviral plasmids   |
| Components:               | ANGPTL4 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 51129). 5µg purified plasmid DNA per construct<br>29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.   |
| RefSeq:                   | <a href="#">NM_001039667</a> , <a href="#">NM_016109</a> , <a href="#">NM_139314</a> , <a href="#">NR_104213</a> , <a href="#">NM_139314.1</a> , <a href="#">NM_139314.2</a> , <a href="#">NM_001039667.1</a> , <a href="#">NM_001039667.2</a> , <a href="#">BC023647</a> , <a href="#">BC023647.2</a> , <a href="#">NM_001039667.3</a> , <a href="#">NM_139314.3</a>   |
| UniProt ID:               | <a href="#">Q9BY76</a>  |
| Summary:                  | This gene encodes a glycosylated, secreted protein containing a C-terminal fibrinogen domain. The encoded protein is induced by peroxisome proliferation activators and functions as a serum hormone that regulates glucose homeostasis, lipid metabolism, and insulin sensitivity. This protein can also act as an apoptosis survival factor for vascular endothelial cells and can prevent metastasis by inhibiting vascular growth and tumor cell invasion. The C-terminal domain may be proteolytically-cleaved from the full-length secreted protein. Decreased expression of this gene has been associated with type 2 diabetes. Alternative splicing results in multiple transcript variants. This gene was previously referred to as ANGPTL2 but has been renamed ANGPTL4. [provided by RefSeq, Sep 2013] |
| shRNA Design:             | These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .  |


[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).