

## Product datasheet for **TL318038**

### **C19orf46 (SYNE4) Human shRNA Plasmid Kit (Locus ID 163183)**

#### **Product data:**

<b>Product Type:</b>	shRNA Plasmids
<b>Product Name:</b>	C19orf46 (SYNE4) Human shRNA Plasmid Kit (Locus ID 163183)
<b>Locus ID:</b>	163183
<b>Synonyms:</b>	C19orf46; DFNB76; KASH4; Nesp4
<b>Vector:</b>	pGFP-C-shLenti (TR30023)
<b>E. coli Selection:</b>	Chloramphenicol (34 ug/ml)
<b>Mammalian Cell Selection:</b>	Puromycin
<b>Format:</b>	Lentiviral plasmids
<b>Components:</b>	SYNE4 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 163183). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
<b>RefSeq:</b>	<a href="#">BC052573</a> , <a href="#">NM_001039876</a> , <a href="#">NM_001297735</a> , <a href="#">NM_153233</a> , <a href="#">NM_001039876.1</a> , <a href="#">NM_001039876.2</a> , <a href="#">NM_001297735.1</a> , <a href="#">NM_001297735.2</a> , <a href="#">BC052573.1</a> , <a href="#">BC038360</a> , <a href="#">NM_001039876.3</a> , <a href="#">NM_001297735.3</a>
<b>UniProt ID:</b>	<a href="#">Q8N205</a>
<b>Summary:</b>	This gene is a member of the nesprin family of genes, that encode KASH (Klarsicht, Anc-1, Syne Homology) domain-containing proteins. In addition to the KASH domain, this protein also contains a coiled-coil and leucine zipper region, a spectrin repeat, and a kinesin-1 binding region. This protein localizes to the outer nuclear membrane, and is part of the linker of nucleoskeleton and cytoskeleton (LINC) complex in the nuclear envelope. LINC complexes are formed by SUN (Sad1, UNC-84)-KASH pairs, and are thought to mechanically couple nuclear components to the cytoskeleton. Mutations in this gene have been associated with progressive high-frequency hearing loss. The absence of this protein in mice also caused hearing loss, and changes in hair cell morphology in the ears. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Aug 2015]
<b>shRNA Design:</b>	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact <a href="mailto:techsupport@origene.com">techsupport@origene.com</a> . If you need a special design or shRNA sequence, please utilize our <a href="#">custom shRNA service</a> .



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).