

Product datasheet for **TL318063**

DNAH10 Human shRNA Plasmid Kit (Locus ID 196385)

Product data:

Product Type:	shRNA Plasmids
Product Name:	DNAH10 Human shRNA Plasmid Kit (Locus ID 196385)
Locus ID:	196385
Synonyms:	DNAH10 variant protein; dynein, axonemal, heavy chain 10; dynein, axonemal, heavy polypeptide 10; FLJ38262; FLJ38262, FLJ43486, FLJ43808, KIAA2017; FLJ43486; FLJ43808; KIAA2017
Vector:	pGFP-C-shLenti (TR30023)
E. coli Selection:	Chloramphenicol (34 ug/ml)
Mammalian Cell Selection:	Puromycin
Format:	Lentiviral plasmids
Components:	DNAH10 - Human, 4 unique 29mer shRNA constructs in lentiviral GFP vector(Gene ID = 196385). 5µg purified plasmid DNA per construct 29-mer scrambled shRNA cassette in pGFP-C-shLenti Vector, TR30021, included for free.
RefSeq:	NM_001083900 , NM_207437 , NM_207438 , NM_207437.1 , NM_207437.2 , NM_207437.3 , NM_207438.1 , BC144575 , BC150622
UniProt ID:	Q8IVF4
Summary:	Dyneins are microtubule-associated motor protein complexes composed of several heavy, light, and intermediate chains. The axonemal dyneins, found in cilia and flagella, are components of the outer and inner dynein arms attached to the peripheral microtubule doublets. DNAH10 is an inner arm dynein heavy chain (Maiti et al., 2000 [PubMed 11175280]). [supplied by OMIM, Mar 2008]
shRNA Design:	These shRNA constructs were designed against multiple splice variants at this gene locus. To be certain that your variant of interest is targeted, please contact techsupport@origene.com . If you need a special design or shRNA sequence, please utilize our custom shRNA service .


[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that the sequences in the shRNA expression cassettes are verified to correspond to the target gene with 100% identity. One of the four constructs at minimum are guaranteed to produce 70% or more gene expression knock-down provided a minimum transfection efficiency of 80% is achieved. Western Blot data is recommended over qPCR to evaluate the silencing effect of the shRNA constructs 72 hrs post transfection. To properly assess knockdown, the gene expression level from the included scramble control vector must be used in comparison with the target-specific shRNA transfected samples.

For non-conforming shRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the shRNA kit. To arrange for a free replacement with newly designed constructs, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled shRNA control (Western Blot data preferred).